

LipoFreeze

***A Modern
Approach to
Cryolipolysis***

by



Innovating for Life



Welcome to the Future of Non-Surgical Fat Reduction

For years, the battle against stubborn fat required invasive procedures and lengthy recovery times. Today, science has delivered a revolutionary solution: **Cryolipolysis**. This groundbreaking, non-invasive therapy uses controlled cooling to target and eliminate fat cells permanently, without surgery or downtime. It's the safe, effective, and proven answer for achieving the contoured physique you desire. Welcome to a new era of aesthetic medicine.

Your Journey to a Sculpted Practice Starts Here




The premier **Cryolipolysis device** engineered for superior performance and remarkable results. Designed and manufactured by **"Health Technology"**

A pioneer in Egyptian medical equipment, **LipoFreeze** combines global technological standards with deep local expertise. We understand the unique needs of the regional market and have built a device that isn't just effective, but also accessible, reliable, and backed by unparalleled support.


LipoFreeze is more than a device; it's your trusted partner in building a successful body contouring practice.




The Science of Slimming: How LipoFreeze Works




Before



Cold Exposure



Crystallization



After

The principle is elegant yet powerful: fat cells (adipocytes) are more susceptible to cold than surrounding tissues. The LipoFreeze applicator delivers precisely controlled cooling to the targeted fat layer.

Crystallization: The cold exposure causes the fat cells to undergo natural apoptosis (programmed cell death).

Natural Elimination: Over the following weeks, your body's lymphatic system naturally processes and permanently removes these dead fat cells.

Visible Contouring: The result is a gradual, natural-looking reduction in the fat layer, with optimal results visible within 2 to 4 months.

Why LipoFreeze Stands Apart: Our Unbeatable Advantages

While others compete on price alone, we build on a foundation of quality and commitment. Our Unique Selling Propositions are the pillars of your success and your clients' satisfaction.

High-Quality Components & Precision Engineering

Every **LipoFreeze** device is built with meticulously sourced, premium-grade components. Our rigorous manufacturing standards ensure optimal cooling performance, consistent results, and exceptional device longevity. We invest in quality so you can invest in your business with confidence.

Competitive Price: Unmatched Value


By leveraging our local manufacturing capabilities in Egypt, we deliver a world-class cryolipolysis device at a highly competitive price point. This strategic advantage makes the LipoFreeze technology accessible to a wider range of clinics and med-spas, maximizing your return on investment without compromising on quality.

Professional After-Sales & Maintenance Services

Your success is our priority. Our commitment begins after your purchase, with a comprehensive after-sales program. This includes professional installation, on-site training for your staff, and routine maintenance checks performed by our certified technicians. We ensure your LipoFreeze device operates at peak performance, always.

Guaranteed Spare Parts Availability

Downtime is lost revenue. That's why we guarantee the immediate availability of all genuine spare parts. Our extensive local inventory means faster service, minimal interruption to your operations, and the peace of mind that comes from working with a reliable and established manufacturer.



Key Benefits at a Glance

Non-Surgical & Painless: No scalpels, no sutures, no downtime.

FDA-Cleared Technology: A safe and clinically proven method for fat reduction.

Targets Stubborn Fat: Effectively reduces fat pockets resistant to diet and exercise.

Long-Lasting Results: Treated fat cells are gone for good.

Ideal Treatment Areas

LipoFreeze is perfect for:
Abdomen & Waistline
Love Handles & Flanks
Inner & Outer Thighs
Back & Bra Roll Fat
Upper Arms





+2 010 23448873



Health-Technology-Egypt



/HealthTechnologyCo



/co.healthtechnology



/HealthTechnology-medsup



Technical Specifications

| | |
|---------------------------|------------------------------|
| Display: | 8 Inch Color Touch Display |
| Operating Systems: | Host |
| Dimensions: | 44 (W) x 47 (D) x 120 (H) cm |
| Power Supply: | 220 V - 50 Hz |
| Temperature: | Down to -10°c |
| Vaccum: | 4 Steps |
| Treatment Time: | Up to 60 min |
| Heads: | 2 Independant Hand Pieces |



The Health Technology Promise

Choosing **LipoFreeze** isn't just about purchasing a device; it's about investing in a complete solution built on quality, support, and value. You gain a competitive edge with a reliable technology that your clients will trust, backed by a team dedicated to your growth.

In **Health Technology**, we are dedicated to innovation, quality, and comprehensive customer care, All are built into every product we create. *We stand behind our technology, so you can stand behind your treatment.*

Ready to enhance your practice with

LipoFreeze

Contact Us Now to see how you can transform your service offerings and boost your profitability.

www.health-technology.net

Main Office: Address: 18 El Shekh Ali Mahmoud – Mahkama Sqr. – Heliopolis, Cairo Egypt - Tel.: +202 264 341 80

Factory: No. 317 Green Land Complex – Obour City, Egypt - Tel.: +202 448 127 08

Email: health@health-technology.net

