

The SHOCKPRO Family



Innovating for Life





What is Shock Wave Therapy Device?

Shock Wave Therapy Devices are advanced medical equipment designed to deliver controlled acoustic pulses to treat musculoskeletal conditions and chronic pain. These non-invasive devices utilize focused or radial shock waves to stimulate tissue regeneration, break down calcifications, and improve blood circulation in affected areas. Commonly used for conditions like **plantar fasciitis, tendonitis, shoulder calcifications, and sports injuries**, shock wave therapy accelerates healing while minimizing patient downtime. Modern devices feature **adjustable energy settings, ergonomic designs, and smart safety protocols** to ensure precise, effective treatments. Ideal for **physiotherapy clinics, sports medicine centers, and rehabilitation facilities**, shock wave technology offers a **drug-free, surgery-free alternative** for pain management and recovery. Its proven efficacy and growing adoption make it an essential tool in modern therapeutic practice.



SHOCKPRO is a cutting-edge radial shock wave therapy device designed to deliver **fast, non-invasive pain relief** and accelerate tissue healing. Engineered for **physiotherapists and clinics**, it features **adjustable energy levels (1.0 - 5.0 bar with 0.5 bar transition)**, 2 different treatment heads for **versatile applications**, and **smart safety technology** for precise, comfortable treatments. **With up to 2000 pulses per session**, it stimulates **cellular regeneration**, reduces **inflammation**, and **improves mobility**—all with **zero downtime**. Proudly made in Egypt with very competitive prices.

Transforms Recovery:

-  **Faster**
-  **Deeper**
-  **Pain-Free**






Ensures reliable performance, easy maintenance, spare parts continuous availability and local technical support, making it the top choice for rehabilitation centers in Cairo, Alexandria, and beyond.



High Quality Stainless Steel Tube and Bullet



2 Different Size Heads for Multiple Functions

Applications:

- ⦿ Myofascial trigger points
- ⦿ Other disorders of tendon insertions
- ⦿ Muscle & connective tissue activation
- ⦿ Cellulite treatment (stages I – III)
- ⦿ Muscle & connective tissue tightening
- ⦿ Elimination of lymphatic congestion
- ⦿ Body shaping
- ⦿ Regulation of muscle tone
- ⦿ Medical massage
- ⦿ Erectyle Disfunction



Dimensions:

Width: 40 cms
Depth: 40 cms
Height: 50 cms
Weight: 20 kg



Dimensions:

Width: 38.5 cms
Depth: 38.5 cms
Height: 110 cms
Weight: 36.5 kg





Enjoy Working with:

- ④ 7" Touch Color Display
- ④ 7 Pre-stored Ready-made Programs
- ④ 3 Open Programs for customization
- ④ Up to 2000 Pulse/Session
- ④ Up to 21 Hz Frequency
- ④ Up to 5 Bar Guaranteed Steady Power
- ④ 0.5 Bar/Transition
- ④ 2 Different Sizes Treatment Heads
- ④ Smart Safety Technology
- ④ 1 Year Full Warranty
- ④ Spare Parts Availability
- ④ Available Maintenance All the Time



+2 010 23448873



Health-Technology-Egypt



/HealthTechnologyCo



/co.healthtechnology



/HealthTechnology-medsup

www.health-technology.net

Main Office: Address: 18 El Shekh Ali Mahmoud – Mahkama Sqr. – Heliopolis, Cairo Egypt - Tel.: +202 264 341 80
 Factory: No. 317 Green Land Complex – Obour City, Egypt - Tel.: +202 448 127 08
 Email: health@health-technology.net

