

# ***THE SONIC FAMILY***

Ultrasonic Physiotherapy Solutions



**SONIC**  
*PLUS*

**SONIC-DUO**

by  **Health**  
TECHNOLOGY  
SINCE 2008

*Innovating for Life*





# ***THE POWER OF ULTRASOUND IN MODERN PHYSIOTHERAPY***

Ultrasound therapy is a cornerstone of modern rehabilitation, utilizing high-frequency sound waves to promote deep tissue healing, reduce inflammation, and alleviate pain. These acoustic waves generate micromassage and thermal effects at a cellular level, enhancing blood flow, increasing tissue elasticity, and accelerating the body's natural repair processes. It is a trusted, non-invasive modality for treating a wide range of musculoskeletal conditions, from chronic tendonitis and muscle strains to joint contractures and scar tissue management. The key to maximizing these benefits lies in delivering precise, reliable, and targeted treatment which is exactly where The Sonic Family excels.

## ***THE SONIC FAMILY***

A suite of advanced ultrasonic devices engineered for clarity, control, and clinical excellence.

**The Sonic Family** offers a scalable solution for every practice, from compact portability to comprehensive multi-modality integration, all while providing the industry's best price-to-quality solution.

Built with the practitioner in mind, each device combines intuitive operation with robust performance, offering complete therapeutic versatility with both Continuous and Pulsed operating modes, adjustable in precise %10 increments from %10 to %90. This unparalleled control is complemented in our advanced models by a unique combination of Ultrasonic and EMS therapy within a single, integrated device, maximizing treatment efficacy and practice efficiency.

Engineered to deliver exceptional patient outcomes, **The Sonic Family** is the smartest investment for your clinic's future.





### **1. SONIC PLUS 1 HEAD**

#### **Precision in its Purest Form**

The foundation of our family, the ideal choice for focused therapy on injury sites like shoulder, elbow, or knee.

#### **Key Features:**

**1 MHz Transducer:** for deep tissue (2-5 cm depth).

#### **Treatment Modes:**

- Continuous
- Pulsed 10 - 90%

#### **Treatment Time:**

- 0 - 30 minutes

#### **Power**

#### **Consumption:**

- 40 W / Hour

#### **Weight:**

- 3 Kg

#### **Power Input:**

- 220 V / 60 Hz

#### **Treatment Area:**

- 1 MHz: 4 cm



### **2. SONIC PLUS 2 HEADS**

#### **Dual-Frequency Versatility**

This innovative device features two independent transducers, to choose the perfect frequency for each clinical scenario..

#### **Key Features:**

**1 MHz Transducer:** for deep tissue (2-5 cm depth).

**3 MHz Transducer:** for superficial tissue focus (1-2 cm depth)

#### **Treatment Modes:**

- Continuous
- Pulsed 10 - 90%

#### **Treatment Time:**

- 0 - 30 minutes

#### **Power**

#### **Consumption:**

- 40 W / Hour

#### **Weight:**

- 3.5 Kg

#### **Power Input:**

- 220 V / 60 Hz

#### **Treatment Area:**

- 1 MHz: 4 cm
- 3 MHz: 3 cm



### **3. SONIC DUO 1 HEAD**

#### **Synergistic Therapy, Unmatched Results**

Combines the deep healing power of ultrasound with the muscle re-education and pain-blocking benefits of Electrical Muscle Stimulation (EMS) saving time and amplifying therapeutic effects.

#### **Key Features:**

**1 MHz Transducer:** for deep tissue (2-5 cm depth).

#### **Treatment Modes:**

- Continuous
- Pulsed 10 - 90%

#### **Treatment Time:**

- 0 - 30 minutes

#### **Power**

#### **Consumption:**

- 40 W / Hour

#### **Weight:**

- 3.5 Kg

#### **Power Input:**

- 220 V / 60 Hz

#### **Treatment Area:**

- 1 MHz: 4 cm

#### **Dual-Channel EMS with Triple Current**

**Options:** Faradic, Low TENS, and High TENS.



### **4. SONIC DUO 2 HEADS**

#### **The Ultimate All-in-One Clinical Solution**

Represents the pinnacle of ultrasonic and electrotherapy technologies. It offers unparalleled flexibility and power with dual-frequency ultrasound and dual-channel EMS.

#### **Key Features:**

**1 MHz Transducer:** for deep tissue (2-5 cm depth).

**3 MHz Transducer:** for superficial tissue focus (1-2 cm depth)

#### **Treatment Modes:**

- Continuous
- Pulsed 10 - 90%

#### **Treatment Time:**

- 0 - 30 minutes

#### **Power**

#### **Consumption:**

- 40 W / Hour

#### **Weight:**

- 4 Kg

#### **Power Input:**

- 220 V / 60 Hz

#### **Treatment Area:**

- 1 MHz: 4 cm
- 3 MHz: 3 cm

#### **Dual-Channel EMS with Triple Current**

**Options:** Faradic, Low TENS, and High TENS.





+2 010 23448873



Health-Technology-Egypt



/HealthTechnologyCo



/co.healthtechnology



/HealthTechnology-medsup



## Why Choose *THE SONIC FAMILY?*

- **Clinical Excellence:** Engineered for effective and reliable therapeutic outcomes.
- **Practitioner-Focused Design:** Durable with clear LCD Display, User Friendly Interface and suitable light weight.
- **High Quality Material:** State-of-the-art components & long-lasting Stainless Steel transducers
- **Scalable Solutions:** A perfect device for every need & every stage of your practice's growth.
- **Dedicated Support:** Backed by our commitment to customer service, maintenance, training, technical support, & affordable available spare parts
- **Competitive Prices:** Engineered for effective and reliable therapeutic outcomes

### Ready to transform your practice?

Contact us today to find the perfect *SONIC FAMILY DEVICE* for you.

[www.health-technology.net](http://www.health-technology.net)

Main Office: Address: 18 El Shekh Ali Mahmoud – Mahkama Sqr. – Heliopolis, Cairo Egypt - Tel.: +202 264 341 80

Factory: No. 317 Green Land Complex – Obour City, Egypt - Tel.: +202 448 127 08

Email: [health@health-technology.net](mailto:health@health-technology.net)

